**ACIDO
IALURONICO
HYALURONIC ACID**

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**HYALURONIC ACID DETECTION
(One Step Detection)**

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Hyaluronic acid is produced by connective tissue cells (fibroblasts). The main biological activity of the molecule is expressed in the connective tissue where it performs a structural function and in the synovial fluid where it is responsible for normal water retention, fundamental for the lubrication of the joints. In plasma the half-life of the molecule has been estimated at approximately 5-6 minutes due to a removal process carried out by the receptors of the sinusoidal endothelial cells of the liver and by the enzymatic action of hyaluronidase.

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Due to their high biocompatibility, viscoelasticity, plasticity and permeability, hyaluronic acid and its sodium salt are widely used in the cosmetics industry and are also widely used in various fields of medicine (ophthalmology, orthopaedics, gynecology, surgery, cosmetology). The analytical method described in the European Pharmacopoeia (01/2011:1472 in Ph. Eur. 8.6) for the determination of hyaluronic acid content is an indirect method, while the commonly used instrumental biochemical methods are complex, expensive to run and require technical expertise. As a result, we decided to develop a direct, accurate and simple method for determining the concentration of native and soluble cross-linked hyaluronic acid.

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The kit allows to perform accurate routine analyses, in a single step, of the cosmetics, drugs and medical devices and can be used for stability studies.

The analytical protocol involves the quantification of native hyaluronic acid and the soluble fraction of cross-linked hyaluronic acid.

Due to the easy execution mode, the kit not require technical expertise.

The performance characteristics of the method are:

good repeatability, high sensitivity, low susceptibility to matrix interference.

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
Hyaluronic acid reacts in a weakly acidic solution with bisphenylnaphthylmethane generating a complex whose optical density at 680 nm is proportional to the concentration of hyaluronic acid in the range 0.001 - 0.01 %. Fig. 1 shows the calibration curve.

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Typical straight line Standard Curve. Spectrophotometres vary in their design and performance so this Figure should be considered as a guide only.

Bibliografia

References

- [1] L. H. Chen, Y. Jian, H. Q. Luo et al. *Chinese Chemical Letters* 18 (2007) 1099–1102
- [2] L. H. Chen , S. Liu , H. Q. Luo and X. Hu. *Journal of Chemical and Pharmaceutical Research*, 2014, 6(6): 1695-1698
- [3] Z. Rongqing, G. Siyuan, L. Lin, H. Songqing. *Journal of South China University of Technology*, 2001, 29(11): 55-58
- [4] L. Jinglan, Z. Chunping, L. Zhiku, X. Dongyan, C. Jinli. *Clinical Medicine of China*, 2002, 18(2): 126-127
- [5] S. Wei, G. Ruifang, D. Yaqin, J. Kui. *Chinese Journal of Analytical Chemistry*, 2006, 34(9): 1265-1268
- [6] T. Bitter, H.M. Muir, *Anal. Biochem.* 4 (1962) 330.
- [7] J.L. Li, C.P. Zhao, Z.K. Lin, D.Y. Xi, et al. *Clin. Med. China*. 18 (2002) 126.
- [8] K. Ruckmania Saleem, Z. Shaikha et al. *Journal of Pharmaceutical Analysis*, Vol. 3 Issue 5, Oct. 2013, 324-329.