

PREMIER COLORSCAN INSTRUMENTS
SERVING THE INDIAN & GLOBAL
INDUSTRIAL SPACE FOR MORE
THAN 35 YEARS

PCI DISPERTECH 75

“EMPOWERING Your R&D team, high speed,
highly reliable lab scale disperser”

**Premier Colorscan
Instruments Pvt Ltd**

has been a pioneer in the field of computerized Color Matching Systems for over three decades. We have been working closely with Paints, Pigments, chemicals, dyestuff and other industries and have understood the requirement of the industry.

Having been in the field of Dispersion technology for over two decades, we have put our experience in developing a technologically advanced disperser. It is our pleasure to introduce the **PCI DISPERTECH 75** – A powerful tool for lab scale dispersion.



FOCUS ON QUALITY

Premier engineering team has focused on designing a highly reliable and quality product, A product that meets the demands of manufacturing industry. This reliability translates into high confidence from your R&D team and which then translates into highly quality end product. The design is focuses on practicality of operations. With high quality mechanical components coupled with high quality electronics you are assured long running life of product. The product is priced economically and meets the needs of Dispersing, stirring, Emulsifying, Suspending.



TECHNOLOGY

Attention to details and the ergonomic design means easy of operation. The 7" Touch Screen display is easy to operate with quick options at the touch of a finger. There is immense focus on safety with 2 step safe operation as a standard feature. Sensors, detect hold of the mixing vessel and height of the dispersion blade to ensure that blade is not exposed to the operator without being immersed in the vessel.

- ❖ High Quality motors specially manufactured in the USA for Premier Colorscan Instruments Pvt Ltd.
- ❖ Stainless Steel components.
- ❖ Ergonomic design for ease of use
- ❖ 2 step safety - Height and Clamp.

SPECIFICATION:

Power (Kilo-Watts) (KW)	0.75
Power (HP)	1
Supply (V)	230
Phase	1
Frequency (Hz)	50
Speed (RPM)	0 - 15000
Torque (nM)	0.5 nM
Orientation	Vertical
Dimensions	L 22 x W 21 x H 31 inch
Container Capacity	Options available
Weight	36 kgs



Premier Colorscan Instruments Pvt. Ltd.

- 📍 EL-130, Electronic Zone, M.I.D.C., Mahape, Navi Mumbai 400 710 India.
- ☎ +91-9004054444 | +91-9004054433 | +91-9004054425
- ✉ colorlab@premiercolorscan.com
- 🌐 www.premiercolorscan.com