

PAL-PH PLUS AT \$350.00 LIST PRICE



Options:

NIST Calibration Certificate (Add \$300.00)
RE-38048 Battery Case Cover (Clear) (Add \$22.00)
RE-39409 Carrying Case (Add \$96.00)
RE-39414 JORDAN for PAL (Silicon Cover) (Add \$27.00)
RE-99210 (Standard solution for calibration pH 4.01) (Add \$55.70)
RE-99211 (Standard solution for calibration pH 6.86) (Add \$55.70)
RE-99212 (Standard solution for calibration pH 7.00) (Add \$55.70)
RE-99213 (Standard solution for calibration pH 9.18) (Add \$55.70)
RE-99214 (Standard solution for calibration pH 10.01) (Add \$55.70)

PAL-pH and PAL-Fruits & Jam pH have been merged to become PAL-pH Plus.
This product has been made even more durable and can be used for a longer period of time.

Product information

Model	PAL-pH Plus
Cat.No.	4301
Range	pH0.00~14.00
Resolution	pH0.01
Accuracy	pH±0.10
Calibration	Calibration at 3 points at 4.01, 6.86, 9.18. Calibration at 3 points at 4.01, 7.00, 10.01.
Temperature compensation	10 to 40° C
Sample Volume	0.6mL
Power Supply	Size AAA alkaline battery × 2